

Universal handle for all reloads (SULUs)

Maximum of 25 firings per device

6 rows of staples on all Endo RLC reload codes

The Endo RLC reloads are available in the following sizes:

- 45 and 60mm lengths
- 2.5, 3.5, 4.0 and 4.8 mm staples for each length size

Gold reload (4.0mm staples) available

Up to 45 degrees of articulation in either direction

24mm opening Jaw

Dotted surface reloads

Reduces the risk of site contamination

Improves OR efficiency and reduces inventory codes

Cost effective

Provides superior hemostasis

Handles a wide range of tissue thickness

- The Gold reload can be used for regular/thick tissue and can pass through a 12mm trocar
- Reduces the use of the Green reload (4.8mm staples), which requires a 15mm trocar

Provides more flexibility and improves access to difficult areas

Easier to insert tissue

Increases the clamp force and reduces tissue slippage during firing

Endo RLC Reloads (SULU)



Product Code	Description	Color	Tissue	Open Staple Height (mm)	Closed Staples Height (mm)	Number of Staples	Staple Rows	Qty/Box
ENDO RLC4525 L	45 mm Linear White SULU	White	Vascular/Thin	2.5	1.0	66	6	6
ENDO RLC4525 R	45mm Articulating Whit SULU	White	Vascular/Thin	2.5	1.0	66	6	6
ENDO RLC4535 L	45mm Linear Blue SULU	Blue	Regular	3.5	1.5	66	6	6
ENDO RLC4535 R	45mm Articulating Blue SULU	Blue	Regular	3.5	1.5	66	6	6
ENDO RLC6035 L	60mm Linear Blue SULU	Blue	Regular	3.5	1.5	90	6	6
ENDO RLC6035 R	60mm Articulating Blue SULU	Blue	Regular	3.5	1.5	90	6	6
ENDO RLC6040 L	60mm Linear Gold SULU	Gold	Regular/Thick	4.0	1.8	90	6	6
ENDO RLC6040 R	60mm Articulating Gold SULU	Gold	Regular/Thick	4.0	1.8	90	6	6
ENDO RLC4548 L	45mm Linear Green SULU	Green	Thick	4.8	2.0	66	6	6
ENDO RLC4548 R	45mm Articulating Green SULU	Green	Thick	4.8	2.0	66	6	6

Endo RLC Universal Stapler

Product Code	Description	Shaft Length (mm)	Max Firings	Qty/Box
ENDO RLC S	Universal Short Endoscopic Linear Cutter	65	25	3
ENDO RLC	Universal Standard Endoscopic Linear Cutter	155	25	3
ENDO RLC L	Universal Long Endoscopic Linear Cutter	255	25	3

Cincinnati, USA
 4480 Lake Forest Dr, Suite 414
 Cincinnati, OH 45242, USA
 Tel: +1(513)769-1918 Fax: +1(513)769-1921

Tianjin, China
 4th Floor, Bldg. B, No. 5 of Fourth
 Street, TEDA, Tianjin, 300457, China
 Tel: +86 22-25325428 Fax: +86 22-66200915

Axon Medical Solutions Pvt. Ltd.

INDIA
 2180, Lower Ground Floor, Sector-46
 Gurugram-122002, Haryana (India)
 Ph.: +91-124-2581022, 9599017699, 9311145832
 Email : info@axonmedicalsolutions.com
 Website : www.axonmedicalsolutions.com

<http://www.reachsurgical.com>

Please always refer to the Instructions for Use
 July 26, 2017

