


PRODUCT SPECIFICATION
ENDO III STAPLERS



Product Code	Description	Shaft Length (mm)	Max Firings	QTY/ BOX
REC60A S	Articulating Endo III 60mm - SHORT	305	12	3
REC60A	Articulating Endo III 60mm - STANDARD	355	12	3
REC60A L	Articulating Endo III 60mm - LONG	455	12	3
REC45A S	Articulating Endo III 45mm - SHORT	290	12	3
REC45A	Articulating Endo III 45mm - STANDARD	340	12	3
REC45A L	Articulating Endo III 45mm - LONG	440	12	3
REC60 S	Soild Shaft Endo III 60mm - SHORT	305	12	3
REC60	Soild Shaft Endo III 60mm - STANDARD	355	12	3
REC60 L	Soild Shaft Endo III 60mm - LONG	455	12	3
REC45 S	Soild Shaft Endo III 45mm - SHORT	290	12	3
REC45	Soild Shaft Endo III 45mm - STANDARD	340	12	3
REC45 L	Soild Shaft Endo III 45mm - LONG	440	12	3



REACH ENDO III
Endoscopic Linear Cutting Staplers and Reloads



INDIA Corporate Office:

Office No: 29, 1st Floor,
Sushant Vyapar Kendra,
Sushant Lok-I, Gurugram,
PIN - 122002, Haryana [INDIA]

Tel : +91 95990 17699
+91 93111 45832
+91 73308 13300

☎ : +91 124 4257540

✉: info@axonmedicalsolutions.com
🌐: www.axonmedicalsolutions.com

R&D Center

4480 Lake Forest Dr, Suite 414
Cincinnati
45242 OH
USA

T +(513)769-1918
F +(513)769-1921




Please always refer to the Instructions for Use



PRODUCT SPECIFICATION

ENDO III Reloads (SULU)



Product Code	Description	Color	Tissue	Open Staple Height (mm)	Closed Staple Height (mm)	Number of Staples	Staple Rows	QTY/BOX
REC 60 GRY	60mm Gray SULU	Gray	Vascular/Thin	2.0	0.75	90	6	12
REC 60 WHT	60mm White SULU	White	Vascular/Thin	2.5	1.0	90	6	12
REC 60 BLU	60mm Blue SULU	Blue	Regular	3.5	1.5	90	6	12
REC 60 GLD	60mm Gold SULU	Gold	Regular	3.8	1.75	90	6	12
REC 60 GRN	60mm Green SULU	Green	Regular/Thick	4.1	2.0	90	6	12
REC 60 BLK	60mm Black SULU	Black	Regular/Thick	4.4	2.2	90	6	12
REC 45 GRY	45mm Gray SULU	Gray	Vascular/Thin	2.0	0.75	66	6	12
REC 45 WHT	45mm White SULU	White	Vascular/Thin	2.5	1.0	66	6	12
REC 45 BLU	45mm Blue SULU	Blue	Regular	3.5	1.5	66	6	12
REC 45 GLD	45mm Gold SULU	Gold	Regular	3.8	1.75	66	6	12
REC 45 GRN	45mm Green SULU	Green	Regular/Thick	4.1	2.0	66	6	12
REC 45 BLK	45mm Black SULU	Black	Regular/Thick	4.4	2.2	66	6	12

FEATURES AND BENEFITS REACH ENDO III




Enhanced clamping force

Prevents unwanted movement and maximizes tissue compression, while preventing unnecessary tissue damage.



Single rotation and articulation knob

Allows rotation and articulation from the same knob, providing ease of use.




Unrestricted, freely articuable tip that rotates to a maximum of 45 degrees in each direction

Provides unparalleled control and precision to access challenging anatomy.




Full range of reloads to cover all types of surgery

- 45 and 60mm lengths.
- Gray / White / Blue / Gold / Green / Black reloads for each length.



Lower profile green reloads (4.1mm staples)

Used for handling regular/thick tissue while still able to pass through a 12mm trocar.



Used reload prevention safety

Ensures that used reloads cannot be reloaded and used, safeguarding patients.



New knife with every firing

Ensures clean cuts, prevents cross-contamination and tumor metastasis.

