

REACH[®] PPH 33


Procedure Set for
Hemorrhoid and Prolapse



Product Overview

The hemorrhoid and prolapse stapler is designed to provide treatment for prolapse and hemorrhoids

CE 0197
510K Number: K120179



More Staples, Better Hemostasis Larger Cavity, Better Resection

Features and Benefits

Features	Benefits
Adjustable staple height technology (0.8 mm to 1.5 mm)	Allows the surgeon to control the level of compression applied to tissue
Reach PPH 33 has 32 staples	Better hemostasis, better specimen
Reach PPH 33 has 4 suture thread ports	Easier retrieval of suture through instrument housing
Fixed anvil	Reduces risk of accidental detachment
Reach PPH 33 has larger housing cavity	Accommodates more tissue
The length of central rod is 63 mm	Clearer sight and better view during operation

Specification

Product Code	Description	Staple Rows	Number of Staples	Staple Leg Length (mm)	Closed Staple Height Range (mm)	Cutting Area Diameter (mm)	Qty/Box
PPH 33	Procedure Set for Prolapse and Hemorrhoid	2	32	3.5	0.8-1.5	24.4	3



India Corporate Office

AXON Medical Solutions Pvt.Ltd.

Office No: 29, 1st Floor,
Sushant Vyapar Kendra, Sushant Lok-I,
Gurugram, PIN - 122002, Haryana [INDIA]

☎ : +91 95990 17699 / 93111 45832 / 73308 13300

☎ : +91 124 4257540

✉ : info@axonmedicalsolutions.com

🌐 : www.axonmedicalsolutions.com

R & D Center

4480 Lake Forest Dr, Suite 414
Cincinnati, 45242 OH, USA
T +1(513)769-1918
F +1(513)769-192

