

To be used in abdominal, gynecologic, thoracic and pediatric surgeries for transection, resection and creation of anastomosis.



Features & Benefits :

- » New knife blade with each firing,Ensures smooth, precise and safe cutting performance while reducing the possibility of tumor metastasis.
- » Tissue Gap Control Mechanism,Ensures consistent staple formation from the proximal end to the distal end.
- » Longer staple line than cutting edge, Prevents blood leakage
- » Protect the medical staff from blade cutting; prevent loading or firing used loading unit.
- » Double-side firing knob,Can be fired from either side, more convenient to operate.
- » Shortcut release button, may easily and quickly open the forks.

Product Description :

Product Code	Description	Packaging
RLC6038S	Linear Cutter 60-3.8	3 pcs/case
RLC6048S	Linear Cutter 60-4.8	3 pcs/case
RLS4535S	Linear Cutter 80-3.8	3 pcs/case
RLS4548S	Linear Cutter 80-4.8	3 pcs/case
RLS6035S	Linear Cutter 100-3.8	3 pcs/case
RLS6048S	Linear Cutter 100-4.8	3 pcs/case
RLS3035L	Linear Cutter 60-3.8 SULU	6 pcs/case
RLS3048L	Linear Cutter 60-4.8 SULU	6 pcs/case
RLS4535L	Linear Cutter 80-3.8 SULU	6 pcs/case
RLS4548L	Linear Cutter 80-4.8 SULU	6 pcs/case

RLS6035L	Linear Cutter 100-3.8 SULU	6 pcs/case
RLS6048L	Linear Cutter 100-4.8 SULU	6 pcs/case