

To be used in abdominal, gynecological and thoracic surgical procedures for resection and transection of tissue.



Features & Benefits :

- » Optimized dumbbell anvil pocket provides better staple formation.
- » Strong anvil reduces the risk of deflection of anvil and SULU.
- » Two-stage closure function allows the reposition of the tissue.
- » Cutting-guide of the instrument facilitates the safe and precise transection.
- » Smart handle allows surgeon to identify easily if the stapler has been fired.
- » Longer and slimmer shaft provides better access and visibility, facilitating TME and LAR surgeries with the reduction of cost.
- » Inlaid rubber in the handles to improve better control of the instrument.
- » Automatic or manual pin placement for ease of use.
- » Finger pads design for easy unloading or loading the SULU.
- » 3.5 and 4.8 cartridges are interchangeable, providing 30, 45 and 60 SULU to meet the need of different length of tissue.

Product Description :

| Product Code | Description | Length of Staple Line | Color Coding | Packaging |
|---------------------|----------------------------|------------------------------|---------------------|------------------|
| RLS3035S | Linear Stapler 30-3.5 | 30mm | Blue | 3 pcs/case |
| RLS3048S | Linear Stapler 30-3.5 | 30mm | Green | 3 pcs/case |
| RLS4535S | Linear Stapler 45-3.5 | 45mm | Blue | 3 pcs/case |
| RLS4548S | Linear Stapler 45-4.8 | 45mm | Green | 3 pcs/case |
| RLS6035S | Linear Stapler 60-3.5 | 60mm | Blue | 3 pcs/case |
| RLS6048S | Linear Stapler 60-4.8 | 60mm | Green | 3 pcs/case |
| RLS3035L | Linear Stapler 30-3.5 SULU | 30mm | Blue | 3 pcs/case |
| RLS3048L | Linear Stapler 30-4.8 SULU | 30mm | Green | 3 pcs/case |
| RLS4535L | Linear Stapler 45-3.5 SULU | 45mm | Blue | 3 pcs/case |
| RLS4548L | Linear Stapler 45-4.8 SULU | 45mm | Green | 3 pcs/case |
| RLS6035L | Linear Stapler 60-3.5 SULU | 60mm | Blue | 3 pcs/case |
| RLS6048L | Linear Stapler 60-4.8 SULU | 60mm | Green | 3 pcs/case |